

## Resistance list:

Medium	UNIVERSAL Lift pump LP	UNIVERSAL Lift pump LP
	Order Ref. No. S 25 – 34 625 12 S 200 – 34 625 29	Order Ref. No. C 25 – 34 625 50 C 200 – 34 625 67
Acetone	O	X
Methyl cyanide	O	X
Diacetone alcohol	O	X
Amines	O	X
Ammonium carbonate	O	X
Amyl acetate	O	X
Anti-freeze	X	X
Benzene	O	X
Bleaching agents	O	X
Detergents / components	X	X
Ethylene dichloride	O	X
Ethyl acetate	O	X
Ethylene glycol	X	X
Ethylene oxide	O	X
Glycerine	X	X
Paint thinners	O	X
Methanol	O	X
Methyl acetate	O	X
Methyl acrylate	O	X
Methyl alcohol (Methanol)	O	X
Methyl ethyl ketone	O	X
Methyl isobutyl ketone	O	X
Methyl isopropyl ketone	O	X
Caustic soda solution	O	X
Paraffin oil	X	X
Potassium ferricyanide	O	X
Potassium hydroxide	O	X
Soap solution	X	X
Sodium acetate	O	X
Sodium hyposulfate	O	X
Sodium nitrate	O	X
Solvents	O	X
Tetrahydrofurane	O	X
Toluene	O	X

**X = resistant / O = not resistant**